

Initial Sample Screening

Before the high-resolution analysis of the epitope and paratope can begin, CovalX performs High-Mass MALDI analysis on the antibody, antigen and the intact antibody/antigen complex.

The goals of these analyses are to verify:

- 1) The integrity of both the antibody and the antigen
- 2) The possible aggregation of the antibody
- 3) The possible multimerization of the antigen
- 4) The stoichiometry of the antibody/antigen interaction (monovalency, bivalency).

CovalX unique High Mass MALDI Detection




Figure 1. Direct detection of the unbound antibody and antigen as well as the intact complex is first determined using the CovalX unique High Mass MALDI detection systems.

CovalX offers a unique analytical service for epitope mapping based on its patented technology. Our service allows ultra-fast epitope characterization of monoclonal antibodies targeting any type of antigen (No limit in molecular weight). Based on mass spectrometry, our epitope mapping service is faster than any other methodology and allows epitope characterization with high resolution. Our technology is adapted to any type of epitope (linear or conformational) and the use of mass spectrometry combined with cross-linking gives unique information on the structure of the antibody-antigen interactions, including paratope identification. In addition, our methods allows additional insight such as stoichiometry of the interaction (mono- or bi-valency), aggregation and integrity of the antibody and antigen.

Crosslinking Mass Spec Mapping (XL-MS)

After initial sample screening is completed, a selected cross-linker reagent, which is shown to effectively stabilize the complex (high mass MALDI analysis) is used for the epitope study. A 50:50 mixture of deuterated to undeuterated cross-linker is created to provide an opportune mass tag for detecting the linkers location within the sequences of both the antigen and the antibody.

Multiple Enzymes Provide Highest Coverage

After stabilization, multiple enzymes are utilized for comprehensive proteolysis. Trypsin, Chymotrypsin, ASP-N, Elastase and Thermolyse are each utilized individually to ensure the most thorough detection possible. The peptide coverage from the sum of these proteolysis provides the highest sequence coverage available with the most overlapping possibilities.


Finally, the peptide data is matched to their sequences using separate software packages to ensure proper identification of both the paratope and epitope from one data set.

- Nazabal et al. Anal. Chem., 2006, 78(11), 3562-3570
- Bich C et al. Anal. Biochem., 2008, 375:35-45

Why choose CovalX's Epitope Mapping Services?

- Characterization of conformational or linear epitopes
- Determination of both epitope and paratope
- Lowest cost commercial service available
- Only 150µg of each antigen/antibody required
- Four weeks delivery time



Epitope Mapping for Linear or Conformational Epitopes: 4 weeks



Antibody/Antigen Interactions

No limit in molecular weight for either the antibody or the antigen. The method is adapted to any protein-protein interaction.

Selected Cross-linker (50:50 deuterated)



Direct detection of the Immuno-complexes by High-Mass MALDI ToF




High-Mass MALDI ToF MS

Using MALDI ToF MS equipped with High Mass system, it become possible to rapidly analyze non-covalent interactions in a very broad mass range (0-2000 kDa) with nM sensitivity

AB Sciex 5800 MALDI ToF mass spectrometer equipped with CovalX HM4 High-Mass system

nLC-Orbitrap MS analysis of stabilized immuno-complexes



Purified Ab/Ag cross-linked complex

Multi-Enzymatic Proteolysis (5 different enzymes)




Overlapping cross-linked peptides

nLC-Orbitrap MS/MS analysis



Q-Exactive Orbitrap High-Resolution PMF

High-Resolution Epitope map



Data Analysis using Interaction softwares
Detection of cross-link peptides (Ab/Ag)

